

Electronic Column and Floor Scales

M304641 – Floor Scale with Cable Display**M304641-01** – Column Scale**M305643** – Body Analyzing Scale with Cable Display**M305643-01** – Body Analyzing Column Scale

- Large Display – 29 mm – for easy readability of all results
- Automatic zero setting and function test
- Tare and hold function
- Automatic switch off
- Low battery indicator
- BMI (Body Mass Index) function
- 25 personal memory settings for all displayed results
- RJ11 interface for digital height measuring rod MZ10033
- Power supply: mains adapter and batteries

Extra Feature of Floor Scales with Cable Display M304641 and M305643:


- Remote display with a 2 m cable

Extra Feature of Column Scales M304641-01 and M305643-01:


- Including wheels for easy manoeuvring
- Optional with height measuring rod mechanical (MZ10023) or digital (MZ10033)

Extra Features of Analyzing Scales M305643 and M305643-01:

- Bioelectrical Impedance Analysis (BIA) technology
- Analysis of body fat, body water, muscle mass and bone mass
- Analysis of fat free mass
- Athlete/Normal mode selection

M304641-01 –
Electronic Column Scale**Technical Specifications**

Capacity	250 kg
Graduation	100 g
Dimensions Column Scales	392 x 320 x 1,014 mm
Dimensions Floor Scales	320 x 320 x 85 mm
Weight Column Scales	approx. 8.2 kg
Weight Floor Scales	approx. 7.1 kg
Platform size	320 x 320 mm

M305643-01 + MZ10023 –
Body Analyzing Column
Scale with mechanical
height rodM304641-01+ MZ10033 –
Height measuring results
are automatically displayedM305643 –
Body Analyzing Floor Scale

M304641 – Floor Scale