

DSmaref Co., Ltd.

High technology of pneumatic medical systems since 1986



Digital Pneumatic Tourniquet Systems

High technology of pneumatic medical systems



DTS-2000w systems...

- 2 channels of individual settings & operations
- 2 channels (4 Ports for 2 double or single cuffs use)
- Fast settings by Touch Key operation
- LCD screen : Easy checking and cleaning
- Save function
- Self test of pressure and system check
- Self air pressure sensor calibration
- Safety alarm system (cuff, pressure, battery, time)
- Built-in battery (Max life time: 12 ~72 hours)
- Max pressure: 600mmHg
- Voltage: Free voltage (AC 100 - 240V~ 50/60Hz)
- Power consumption: 50W
- Noise dB: Below 60dB
- Max pressure: 600mmHg
- Time range: 5~180 Min.
- Compressor weight: 5.6kg
- Compressor dimension: 225(W) X 185(D) X 240 (H)mm

Digital Pneumatic Tourniquet cuffs

High technology of pneumatic medical systems

REUSABLE CUFFS



DISPOSABLE CUFFS



Digital Pneumatic Tourniquet cuffs

High technology of pneumatic medical systems

ACCESSORIES

DTC-S02	SINGLE 40 X 7cm CUFF
DTC-S04	SINGLE 52 X 7.5cm CUFF
DTC-S05	SINGLE 61 X 9cm CUFF
DTC-S06	SINGLE 80 X 9cm CUFF
DTC-S07	SINGLE 86 X 10cm CUFF
DTC-S08	SINGLE 107 X 10cm CUFF
DTC-C25	CONE SINGLE 70 X 10cm CUFF
DTC-C26	CONE SINGLE 90 X 12cm CUFF
DTC-D04	DOUBLE 57 X 12cm CUFF (5
DTC-D05	DOUBLE 80 X 15cm CUFF
DTS-2000	MAIN HOSE, 200cm
DTS-2000	CUFF HOSE, 15cm
DTS-2000	STAND, 102cm X 54cm
DTS-2000	STAND CONNECTOR