

# PRODUCT INFORMATION



## International Range Electronic Baby Scales



**M10614**

## Baby-/ Toddler Scale

- Removable baby tray
- Large LCD, 25 mm
- Unit function, kg / lbs selectable
- Tare function
- Memory and recall function
- Difference value between 2 weighings
- Power supply: Battery

Helps to control baby's feeding

### Technical Specifications

Capacity	20 kg
Graduation	10 g
Dimensions	560 x 105 x 370 mm
Weight	approx. 3.5 kg
Baby tray	560 x 105 x 275 mm