

TECHNICAL DATA & ORDERING INFORMATION: KANMED° WARMING CABINET						
LARGE CABINET				SMALL CABINET		
ART. NO.:	GE-I380 (ADJUSTABLE FEET) GE-I380-W (WITH WHEELS) GE-I380/50 (50° C) GE-I380-W/50 (50° C)			GE-I380-90CM (ADJUSTABLE FEET) GE I380-90CM-W (WITH WHEELS) GE-I380-90CM/50 (50° C) GE-I380-90CM-W/50 (50° C)		
OUTER SIZE:	HEIGHT 175 cm	WIDTH 66 cm	DEPTH 60 cm	HEIGHT 90 cm	W IDTH 66 cm	DEPTH 60 cm
INNER SIZE:	HEIGHT 140 cm	WIDTH 56 cm	DEPTH 53 cm	HEIGHT 65 cm	W IDTH 56 cm	DEPTH 53 cm
WEIGHT:	About 180 Kg			About 110 Kg		
VOLTAGE:	220–240 V AC /50–60 Hz			220–240 V AC /50–60 Hz		
POWER MAX:	900 Watt			900 Watt		
POWER CONSUMPTION:	Average about 150-200 Watt			Average about 150-200 Watt		
FUSES:	T6,3 A L 250 V			T6,3 A L 250 V		
MAX LOAD PER SHELF/DRAWER:	8/20 Kg			8/20 Kg		
BASKET SIZE:	HEIGHT 10 cm	WIDTH 52 cm	DEPTH 50 cm	HEIGHT 10 cm	W IDTH 52 cm	DEPTH 50 cm
SHELF SIZE:	HEIGHT 2 cm	WIDTH 52 cm	DEPTH 50 cm	HEIGHT 2 cm	W IDTH 52 cm	DEPTH 50 cm
STORING HEIGHT BETWEEN SHELVES:	12 cm					
WARMING CAPACITY:	Will warm 10 pieces GE-4550I5 from 22° C to 40° C in less than 3 hours			Will warm 5 pieces GE-4550I5 from 22° C to 40° C in less than 3 hours		
PRE-SET TEMPERATURE:	42° C ± 1° C or 50° C ± 1° C			42° C ± 1° C or 50° C ± 1° C		
TEMPERATURE RANGE:	35-42° C or 35-50° C			35-42° C or 35-50° C		
VOLUME:	LARGE CABINET 415 L			SMALL CABINET 193 L		
OTHER FEATURES:	Made of stainless steel with no flammable insulation Order with adjustable feet or wheels Prepared for mounting of 10 shelves or drawers Practical usable volume per shelf/basket is 25 L			Made of stainless steel with no flammable insulation Order with adjustable feet or wheels Prepared for mounting of 5 shelves or drawers Practical usable volume per shelf/basket is 25 L		

TECHNICAL DATA & ORDERING INFORMATION: KANMED° OPERATHERM	
CONTROL UNIT	
ART. NO.:	OP3-022 (230 V AC) OP3-115 (115 V AC) OP3-100 (100 V AC)
TEMPERATURE RANGE:	33° C - 39° C
HEATING PADS	
ART.NO./DIMENSION/WEIGHT:	OP3-I04, 104 X 45 CM, 3 KG OP3-I50, 150 X 45 CM, 3,5 KG OP3-063, 63 X 43 CM, 1,5 KG OP3-050, 30 X 50 CM, 1,3 KG

TECHNICAL DATA & ORDERING INFORMATION: KANMED° GEL PADS	
ART. NO./DIMENSION/WEIGHT:	GE-4550I5, 45 X 50 CM, 3,5 KG GE-4363I5, 43 X 63 CM, 4 KG GE-3050I5, 30 X 50 CM, 3 KG

For further information please see www.kanmed.se

KANMED° WARMINGCABINET



A HIGH-QUALITY UNIVERSAL WARMING CABINET THAT HELPS PREVENT PERIOPERATIVE HYPOTHERMIA

PRODUCT: **KANMED° WARMINGCABINET**

Made in Sweden. Stainless steel and high quality components ensure a long and trouble free lifetime.

Universal Warming Cabinet for rapid prewarming of:

- Kanmed Anti Pressure Gel Pads
- Gel Positioning Equipment
- Blankets, Sheets & Towels
- Intravenous & Irrigation fluids

Shelves and Drawers glide easily on sturdy roller bearings making access to items easy.

The Cabinet comes with either adjustable feet or lockable wheels.

Double insulating safety glass makes it easy to see the stored items.

The well insulated door and cabinet makes it very silent and the power consumption is low*.



**Compared to less insulated cabinets you will save more than Euro 200 per year in power consumption!*

The KANMED° WARMINGCABINET is designed to assist in the reduction of patient hypothermia. The Cabinet is mainly used in the operating theatre area for warming intra venous, dialysis & irrigation fluids, Kanmed Gel Pads and other gel positioning equipment.

The Cabinet is also used in the recovery area to warm blankets and other patient bedding material.

The Cabinet will actively assist in maintaining normothermia in the perioperative phase.

The universal KANMED° WARMINGCABINET is made in Sweden using stainless steel and high quality components making the construction very sturdy. A long and trouble free life time is ensured.

The KANMED° WARMINGCABINET warms by circulating warm air. Default maximum temperature is 42°C, but a 50°C Cabinet is also available. The temperature is easy to change between 35°C to 42°C (or 35°C to 50°C).

The KANMED° WARMINGCABINET comes in two sizes; height 175 cm or 90 cm. You can equip it with shelves or drawers, lockable wheels or adjustable feet. The large version accommodates up to 10 shelves or drawers, the small version 5 shelves or drawers.



*Kanmed Warming Cabinet Large
Art. No. GE-1380*



*Kanmed Warming Cabinet Small
Art. No. GE-1380-90CM*

A GREAT COMBINATION

The patient warmer KANMED° OPERATHERM is mostly used in the OR area but can also be practical for warming in cardio vascular x-ray procedures, neonatal areas etc. The Kanmed Operatherm together with Kanmed Anti Pressure Gel Pads, that are prewarmed in the KANMED° WARMINGCABINET, will help you prevent patient hypothermia efficiently.



KANMED° OPERATHERM



KANMED° GELPAD



KANMED° WARMINGCABINET