



MACHAN INTERNATIONAL CO., LTD.

No. 59-7, Shan-jiao Lane, Shui-mei Village,

Waipu, Taichung 43858, Taiwan

Tel: +886-4-2688-7171

Fax: +886-4-2688-5050

sales@bailida-medical.com

www.bailida-medical.com

CE ISO 9001:2000



MI SERIES

MEDICAL FURNITURE & HOSPITAL EQUIPMENT

MINI CARTS

The mini cart is designed with smaller overall dimensions ideal for areas with limited space. It provides convenient storage and is easily moved anywhere. Central key lock control also provides the necessary security on this lightweight and versatile unit.



Features

- ABS top & bumpers
- Central key lock to secure all drawers



- A keyboard drawer is available




- Magnetic label holders indicate drawer contents conveniently

Powder paint

4 main colours, other colour options available

- Anti-static paint available
- Anti-static and anti-bacterial paints available

4 main colours



Other available colours



Drawer specifications



3" Drawer
W 399 x D 398 x H 76 mm (3")



6" Drawer
W 399 x D 398 x H 155 mm (6")



9" Drawer
W 399 x D 398 x H 234 mm (9")

SPECIFICATIONS	MI3105		MI3405		MI3706	
	mm	inch	mm	inch	mm	inch
Drawer 3"-6"-9"	3-1-1		2-2-1		3-2-1	
W	520	20.4	520	20.4	520	20.4
D	540	21.2	540	21.2	540	21.2
H	867	34.1	946	37.2	1025	40.3
H (w/o 4" Casters)	732	28.8	811	31.9	890	35.0