



















Fixed divider set			
 P-FD3	3" fixed divider set in drawer (3 x short + 3 x fixed + 2 x 3" strip)	 P-FT3	3" fixed divider set (3 x short + 3 x fixed)
 P-FD6	6" fixed divider set in drawer (2 x short + 2 x fixed + 2 x 6" strip)	 P-FT6	6" fixed divider set (2 x short + 2 x fixed)
Adjustable divider set (Length adjustable)			
 P-AD3	3" adjustable divider set in drawer (3 x short + 6 x adjustable + 2 x 3" strip + 6 x T clips)	 P-AT3	3" adjustable divider set (3 x short + 6 x adjustable + 6 x T clips)
 P-AD6	6" adjustable divider set in drawer (3 x short + 4 x adjustable + 2 x 6" strip + 4 x T clips)	 P-AT6	6" adjustable divider set (3 x short + 4 x adjustable + 4 x T clips)
3" Individual divider			
 X-FD3 3" Fixed divider 542 x 51 x 4.6mm	 X-SD3 3" Short divider 368 x 51 x 4.6mm	 X-AD3 3" Adjustable divider 237 x 51 x 5.3mm	 X-LH Label holder 40 x 60 x 18mm 20 pcs/pack
6" Individual divider			
 X-FD6 6" Fixed divider 542 x 130 x 4.6mm	 X-SD6 6" Short divider 368 x 130 x 4.6mm	 X-AD6 6" Adjustable divider 237 x 130 x 5.3mm	 XTC T clip 40 x 36 x 17.5mm 10 pcs/pack

March 2011 copyright by BAILIDA



MX SERIES

















PEDIATRIC EMERGENCY CART



MACHAN INTERNATIONAL CO., LTD.

Tel: +886-4-2688-7171 Fax: +886-4-2688-5050
 sales@bailida-medical.com www.bailida-medical.com
 No. 59-7, Shan-jiao Lane, Shui-mei Village, Waipu, Taichung 43858, Taiwan



Fixed divider set			
 P-FD3	3" fixed divider set in drawer (3 x short + 3 x fixed + 2 x 3" strip)	 P-FT3	3" fixed divider set (3 x short + 3 x fixed)
 P-FD6	6" fixed divider set in drawer (2 x short + 2 x fixed + 2 x 6" strip)	 P-FT6	6" fixed divider set (2 x short + 2 x fixed)
Adjustable divider set (Length adjustable)			
 P-AD3	3" adjustable divider set in drawer (3 x short + 6 x adjustable + 2 x 3" strip + 6 x T clips)	 P-AT3	3" adjustable divider set (3 x short + 6 x adjustable + 6 x T clips)
 P-AD6	6" adjustable divider set in drawer (3 x short + 4 x adjustable + 2 x 6" strip + 4 x T clips)	 P-AT6	6" adjustable divider set (3 x short + 4 x adjustable + 4 x T clips)
3" Individual divider			
 X-FD3 3" Fixed divider 542 x 51 x 4.6mm	 X-SD3 3" Short divider 368 x 51 x 4.6mm	 X-AD3 3" Adjustable divider 237 x 51 x 5.3mm	 X-LH Label holder 40 x 60 x 18mm 20 pcs/pack
6" Individual divider			
 X-FD6 6" Fixed divider 542 x 130 x 4.6mm	 X-SD6 6" Short divider 368 x 130 x 4.6mm	 X-AD6 6" Adjustable divider 237 x 130 x 5.3mm	 XTC T clip 40 x 36 x 17.5mm 10 pcs/pack

March 2011 copyright by BAILIDA



MX SERIES

PEDIATRIC EMERGENCY CART



MACHAN INTERNATIONAL CO., LTD.


Tel: +886-4-2688-7171 Fax: +886-4-2688-5050

sales@bailida-medical.com www.bailida-medical.com

No. 59-7, Shan-jiao Lane, Shui-mei Village, Waipu, Taichung 43858, Taiwan

MX SERIES MX37EMGP

PEDIATRIC EMERGENCY CART



- Cream >
- Pink/Red >
- Purple >
- Yellow >
- White >
- Light Blue >
- Orange >
- Green >
- Cream >



MX37EMGP2

- Metal with powder coated
- ABS top, bumper & sliding table
- Aluminum frames to protect the cart body
- Drawer : 3" x 8 + 6" x 1, central key lock
- Drawers are Color Coded using for emergency pediatric care
- Cart dimension : W690 x D518 x H1085 mm (w/5" casters)
- Various accessories and drawer divider sets are available
- Caster : 5" (125mm) Iron casters (standard), (2 w/brake, 1 free, 1 w/direction control)
- Lock options : One locking bar for all drawers (MX37EMGPI as the cover)

Individual for each drawer (MX37EMGP2)



ABS top surface
Light, easy to clean.



Sliding table top extension
Easy to use, space saving,
max load 3kg



Push handle
User friendly design,
height adjustable.



Central key lock
secures all drawers

Caster options

Size: 5" (ø125mm)



Twin wheel
TPR caster



Plastic medical
TPR caster



Dustproof medical
TPR caster



Iron medical
TPR caster (standard)

ACCESSORY LIST


Top accessories



VIV
IV pole




VDS
Defibrillator shelf
355 x 355 ~ 415mm




VTB-2
Tilt bin organiser (S)
603 x 106 x 253mm



VTB-1
Tilt bin organiser (L)
603 x 123 x 153mm




P-AB
Accessory bridge
800mm



VBS-I & VPR-I
Bulk storage bins (L) & back rail
75 x 105 x 190mm



VGD & VPR-I
Glove dispenser & back rail
160 x 135 x 95mm




VBH & VPR-I
Bottle holder & back rail
355 x 95 x 80mm

Rear accessories



VCPR & VHK
CPR cardiac board & hooks
440 x 510 x 10mm



VMH
Magnetic label holder
L80 x W50 x T5mm
6pcs/pack



VPT
Pull-tite seal in pink
w/ serial no., 50pcs/pack



VPS
Plastic seal in red
w/ serial no., 50pcs/pack

Others

Side accessories



VGD & VPR-2
Glove dispenser & side rail



VSC & VPR-2
Sharps container & side rail



VBH & VPR-2
Bottle holder & side rail



VMS & VPR-2
Multi-storage (L) & side rail



VWC & VPR-2
Waste container w/ lid & side rail
255 x 170 x 315mm



VWB & VPR-2
Waste bin & side rail
265 x 255 x 460mm



VGB & VPR-2
Garbage bin & side rail
255 x 170 x 315mm



VFH & VPR-2
File holder & side rail
571 x 397 x 150mm



VO2 & VPR-2
Oxygen tank holder & side rail
145 x 125 x 403mm



VCH & VPR-2
Catheter holder & side rail
100 x 100 x 480mm



VSH & VPR-2
Suction machine holder & side rail
213 x 395 x 115mm



VDH & VPR-2
Documents holder & side rail
330 x 160 x 255mm