

# Connecting

# EDAN

from

heart

to

heart



SonoTrax Series Ultrasonic Pocket Doppler's superior sensitivity, high-fidelity sound, recording option and optional rechargeable battery ensure the doppler family to meet the routine examination requirements of obstetricians and midwives.

## SonoTrax

### Ultrasonic Pocket Doppler

#### Features

- Easy-to-use
- Compact design
- Interchangeable waterproof probe
- Built-in audio recorder
- Large backlight LCD display
- Accurate fetal heart detection with high-fidelity sound
- Support 2MHz, 3MHz, 4MHz, 5MHz, 8MHz waterproof probes
- Automatic and manual counting of FHR
- Rechargeable battery optional
- Automatic power-off

# SonoTrax

## Ultrasonic Pocket Doppler



Compact Li-ion Battery Charger



Backlight LCD



Built-in Recorder



Interchangeable Waterproof Probe



Rolling Stand



### Specifications

**Display:** 45x25mm LCD display

**Heart Rate Range:** 50~210bpm

**Resolution:** 1bpm

**Recording Time:** 240s

**Audio Output Power:** 0.5W

**Auto Power-off:** 1 minute with no signal or operation

**Size:** 32 (Depth) x 85 (Width) x 138 (Height) mm

**Weight:** 290±5g (including one battery)

Model	Interchangeable Probes	FHR Display	Built-in Audio Recorder	Backlight	Waterproof Probe	Probe Auto-detection Display	Battery Status Display	Low Battery Indication	AC Recharge
<i>SonoTrax Lite</i>	Yes	No	No	No	Yes	No	No	Yes	No
<i>SonoTrax Basic</i>	Yes	Yes	No	No	Yes	Yes	No	Yes	No
<i>SonoTrax Basic A</i>	Yes	Yes	No	Yes	Yes	Yes	Yes	Yes	No
<i>SonoTrax Pro</i>	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No
<i>SonoTrax II</i>	Yes	Yes	No	Yes	Yes	Yes	Yes	Yes	Yes
<i>SonoTrax II Pro</i>	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
<b>Optional accessories</b>	Rechargeable Li-ion battery, charger, earphone, carry bag, 9V alkaline battery, adapter, Ni-MH battery, probe (2MHz, 3MHz, 4MHz, 5MHz, 8MHz)								

# EDAN

Edan Instruments, Inc. www.edan.com.cn E-mail: info@edan.com.cn  
 Edan Headquarter: 3/F - B Nanshan Medical Equipments Park, Nanshai Rd 1019#, shekou, Nanshan  
 Shenzhen, 518067 P.R. China Tel +86-755-26898326 Fax +86-755-26898330



All rights reserved.

Features and specifications are subject to change without notice

ENG-CTG-SonoTrax-V1.5-20100618