

**Optional accessories available for all models of Victoria:**

Foot switch*	11-5102
Change-over valve	11-5107
Suction jar holder for 2 l bottle for eurorail	11-5108
Suction jar holder for 2 l bottle for side of the unit	11-5113
Suction jar holder for 4 l bottle for eurorail	11-5109
Suction jar holder for 4 l bottle – double	11-5114
Suction circuit holder	11-5110
Storage basket	11-5105
Silicone suction tubing with suction interrupter (1,5 m)	305-888-003
Silicone tubing 8 × 14 mm (25 m)	232-008-014
PVC tubing 8 × 12 mm (20 m)	231-008-003
Holder for catheter container, double	11-5112
Holder for catheter container, single	11-5111
Catheter container without holder	111-888-021
Limiting water valve	11-5106
Suction jar 2 l, screw cover – polysulfone	111-888-288
Suction jar 2 l, press-on cover – polycarbonate	111-888-201
Suction jar 2 l, screw cover – polycarbonate	111-888-200
Suction jar 4 l, press-on cover – polycarbonate	111-888-400
Disposable secreta bag MONOKIT®	114-087-001
Stop valve	111-008-001
Microbiological suction filter MSF(20 pcs/pack)	271-022-001
Integrated catheter holder	11-5111
Mobile stand	11-5104
Mobile stand with large rail	11-5119
Mobile stand with rail 30 × 10	11-5123
Suction jar holder 2 l bottle for rail 30 × 10	11-5124
FLOVAC 2 l jar	970-010-212
FLOVAC 2 l bag	000-036-011
Shelf	11-5116
Eurorail - Victoria Portable	11-5127

\*)Additional installation is only possible in the works.

## Suction unit Victoria



Disposable secreta bag MONOKIT



0,5 - 4 l polycarbonate suction jars autoclavable



New eurorail for placing a large number of optional accessories



MSF microbiological suction filter (20pcs/box)



Easy and safe handling



Change-over valve



IEC power outlet. It is possible to order unit with IEC power outlet.



Air-operated on/off footswitch



**CHEIRÓN a.s.**, Ulrychova 13, 162 00 Prague 6, Czech Republic  
 Branch office: Blatenská 27a, 326 00 Plzeň, Czech Republic  
 Tel.: +420 377 590 433, 422, Fax: +420 377 590 435  
 E-mail: export@cheiron.eu, www.cheiron.eu

In regard of the continuous research and development Cheirón reserves the right to make changes which may affect the information contained in this brochure without prior notice.

CH-06-04-091208\_4003

**Makes less noise than your PC!**

A special system of noise and vibration dampening for an extremely quiet and smooth operation with an incredibly low noise level

**of 39,4 dB.**

# Suction unit Victoria

# Cheirón®

[www.cheiron.eu](http://www.cheiron.eu)



Vacuum gauge in kPa, mmHg and Bar

Illuminated on/off switch

Vacuum regulator

Change-over valve shift the vacuum from one bottle to the other with a single move - optional

Quick-Vac for easy, fast and safe hose connection

The rail ensures maximum flexibility you can attach up to four jars of any size and/or any other accessories needed

- Extremely silent run
- Maintenance-free aggregate
- Reliable and vibration-free run
- Designed for continuous run
- Low maintenance requirements
- 6 year warranty on the aggregate

Non-marking casters with safety brakes.  
Antistatic casters - optional



### Victoria Versa

The most popular high performance model of Victoria especially suitable for wards, operating rooms, recovery rooms, ICUs, obstetrics and endoscopic applications, etc.

### Victoria Lipos

High performance solution designed particularly for liposuction and operating rooms.



### Victoria Portable

Quiet and effective suction unit for medical practices, wards and endoscopic applications, as well suitable for fast suction in any situation.

### Victoria Portable with eurorail

Advantage of the unit is that instead of regular holder, unit is equipped with eurorail 30 x 10 mm this allows you to use your holders and jars.



### Victoria Thorax

A suction unit specially designed for common use as well as for thoracic and abdominal surgery drainage and intensive care.

The system can be used for long-term up to 30 days continuous drainage. Thanks to a unique patented solution Victoria Thorax enables full suction output even during drainage.

Type	Portable	Versa	Thorax	Lipos
<b>Order no:</b>	11-1110	11-1120	11-1130	11-1140
Suction jar 2 l with press on cover	•	•		•
Suction jar 2 l with screw cover			• 2 ks	
Suction jar holder for 2 l bottle for the side of the unit	•			
Suction jar holder for 2 l bottle for eurorail		•	• 2 ks	•
PVC suction circuit with stop valve	•	•	•	•
Quick-Vac hose connection	•	•	•	•
Mobile stand		•	•	•
Limiting water valve			•	
Suction capacity nominal +/- 5 %	43 l/min	43 l/min	43 l/min	50 l/min
Suction capacity effective +/- 5 %	38 l/min	38 l/min	38 l/min	45 l/min
Vacuum	0-90 kPa (0-680 mmHg)	0-90 kPa (0-680 mmHg)	0-90 kPa (0-680 mmHg)	0-93 kPa (0-700 mmHg)
Drainage	0 - 40 cmH <sub>2</sub> O	0 - 40 cmH <sub>2</sub> O	0 - 40 cmH <sub>2</sub> O	0 - 40 cmH <sub>2</sub> O
Noisy level	LAeq 39,4 dB	LAeq 39,4 dB	LAeq 39,4 dB	LAeq 39,4 dB
Input	110 VA	110 VA	110 VA	110 VA
Dimensions (w x h x l) in millimeters	440 x 300 x 310	480 x 900 x 460	480 x 900 x 460	480 x 900 x 460
Dimensions (w x h x l) in inches	20,1 x 11,8 x 13,4	18,1 x 18,9 x 39	18,1 x 18,9 x 39	18,1 x 18,9 x 39
Weight (without accessories)	12 kg (26,5 lb)	21 kg (46,3 lb)	24 kg (52,9 lb)	21 kg (46,3 lb)

[www.cheiron.eu](http://www.cheiron.eu)