



YP-90
INFANT INCUBATOR



YP-90B
INFANT INCUBATOR

Brief introduction

- Air mode control, set temperature and air temperature are displayed separately;
- Natural air flow humidity;
- Inclination of bassinet is adjustable;
- Failure alarm indication;
- RS-232 connector
- YP-90B is configured with fixed Neonate Bilirubin Phototherapy Equipment model: XHZ-90 (with blue fluorescent lamps)
- YP-90B is configured with fixed Neonate Bilirubin Phototherapy Equipment model: XHZ-90L (with LED bulbs) as optional

Endless love keeps the new life safe

Standard configuration

YP-90: The hood, the base, the controller, I.V pole, shelf.

YP-90B: Based on YP-90, with Neonate bilirubin phototherapy equipment (fixed model)

Optional parts: > 37°C temperature Override Mode function, Vertical Height Adjustment Stand(VHA Stand).

Specification

- Power supply: AC220V-230V/50Hz
- Max Power Output: 850VA
- Mode of control: Air mode controlled by micro-computer
- Air temperature control range: 25°C-37°C
(Override mode: 37-39°C)
- Temperature variability: $\leq 0.5^{\circ}\text{C}$
- Temperature uniformity: $\leq 0.8^{\circ}\text{C}$ (mattress is on the horizontal position)
- Inclination of the bassinet: $\pm 12^{\circ}$
- Noise inside hood: $\leq 55\text{dB(A)}$ (Environment noise $\leq 45\text{dB(A)}$)
- Failure alarm: over temperature alarm, deviation alarm, sensor failure alarm, fan motor alarm, power failure alarm, system alarm and so on.

Environmental temperature

Operating range: 20 °C~30°C

(Set temperature should be higher 3°C than the environment temperature)

Ambient air movement rate: $\leq 0.3\text{m/s}$

Package

Package: each unit is packed as one carton

Size: 108.5cm x 64cm x 112.5cm

Gross Weight: 80Kg (for YP-90)

90Kg (for YP-90B)



Official South
African Distributor

NINGBO DAVID MEDICAL DEVICE CO.,LTD.

Sales Center: #100 Jinghua Rd, High-tech industrial development zone, Ningbo China

Tel: 0086-574-87800008 87800009

Fax: 0086-574-87803714

P.C.: 315040

Http://www.nbdavid.com

E-mail: sales@chinadavid.cn